

# Isotopic gases

## Technical data sheet

Our standard cylinder holds 1 L at 12 bar and features a 7/16" DIN valve connection. Bigger cylinders and suitable valves are available on request.



### Valve portfolio:

#### Fine Control Valve without Pressure Gauge

- Inlet pressure: max. 15 bar
- Outlet pressure: same as inlet pressure



#### Fine Control Valve with Pressure Gauge

The pressure gauge allows monitoring of the filling level of the pressure canister.

- Inlet pressure: max. 15 bar
- Outlet pressure: same as inlet pressure

